



# UltraSPX-S200/S220

Scanning-Probe X-ray Microscope for Synchrotron Facilities



## Key Benefits

Structural, elemental, chemical and crystallographic information on the nanoscale

Ultra-high spatial resolution to 30 nm

Fluorescence, absorption, XANES, phase and / or crystallographic contrast

Open system architecture offers flexible configuration

Cryogenic sample handling to minimize the effects of radiation damage for organic specimens

## Ultra-high resolution analysis and characterization with a scanning nanoprobe

The UltraSPX takes advantage of the coherent X-ray beams available at modern synchrotron light sources to bring established X-ray analytical techniques, such as X-ray fluorescence, spectroscopy, and diffraction, down to the nanometer scale. Xradia's high-performance X-ray optics produce a focused beam down to 30 nm. This X-ray probe is raster scanned across the sample to produce highest-resolution maps with elemental, chemical, or crystallographic phase information depending on the photon energy and detection modality.

## Application Examples

### • Materials Science

Characterization of nano-materials in areas such as alternative energy, fuel cells, micromechanics, and industrial processes

### • Life Science

Mapping and quantification of trace elements in cells or tissue samples with parts-per-billion sensitivity, providing insight into uptake, function and toxicity in biological systems


### • Geosciences and Environmental

Analysis of specimens in their natural environment, critical for areas such as heavy metal uptake, nuclear waste storage and bioremediation

### • Semiconductor


Probing of properties such as strain or doping concentration that affect functionality of modern ICs on the nanoscale, without destroying the component.

Laser-interferometer technology and precision piezo stages provide sub-5 nm positioning of beam and sample with excellent long-term stability. In the soft X-ray range, typically transmitted intensity is detected as a function of position and photon energy, to produce X-ray Absorption Near-Edge Structure (XANES) image stacks revealing chemical composition and the binding state of elements. At higher X-ray energies, elemental information can be acquired by detecting fluorescence photons with sensitivities as low as parts-per-billion. Additional detection modalities such as a configured transmission detector for phase contrast imaging or a diffraction detector for crystallographic phase mapping can be added easily to the system to provide additional structural information.



Iron Fluorescence of a geo-polymer.

Image courtesy of J. Provis (Univ. Melbourne) and V. Rose (ANL).




Coherent Bragg Diffraction Pattern of a functional nanodevice.

Image courtesy of M. Bedzyk (Northwestern Univ.) and M. Holt (ANL)



Automated Sample Exchange Robot for cryogenic sample handling system.



Diffraction pattern of a thin film sample using a zone plate.

Image courtesy of C. Thompson (NIU), G. B. Stephenson (ANL) and M. Holt (ANL)

# UltraSPX-S200/S220 Scanning-Probe X-ray Microscope for Synchrotron Facilities

## Features

Multiple zone plate optics possible to tune beam flux and resolution to experimental requirements

Integration of KB mirrors and other optics possible


Integrated visible light microscope for easy sample alignment

Cryogenic sample handling with automated sample exchange robot available

System flexibility and open EPICS-based control system for integration of in-situ devices, additional detector and general extensibility

Laser interferometer controlled step- and fly scanning


Sub-5 nm positional stability



## Key Specifications

Systems	UltraSPX-S200	UltraSPX-S220
Energy range	5 to 20 keV	200 to 3000 eV
Best Spatial Resolution	30 nm	30 nm
Sample environment	Air or vacuum	Vacuum
Positional stability	< 5 nm	< 5 nm
Detector options	Fluorescence Photodiode Transmission camera Diffraction camera possible (customer integration)	Phosphor/PMT Photodiode Integration of customer-supplied detectors possible (e.g. pixel detector)
Data acquisition times	Beamline and application dependent	

\* Specifications are typical. Please consult Xradia for custom requirements.



The UltraSPX-S200 / S220 typically requires an undulator beamline. Please consult Xradia for beamline recommendations.

**All specifications subject to change. Please consult Xradia for current specifications.**



**Xradia, Inc.** 4385 Hopyard Road, Pleasanton, CA 94588

T: +1.925.701.3600 F: +1.925.730.4952 sales@xradia.com [www.xradia.com](http://www.xradia.com)

Xradia, VersaXRM, UltraXRM, MicroXCT, Scout-and-Zoom, Scout-and-Scan, and PhaseEnhanced Optics are trademarks of Xradia, Inc. The VersaXRM family of X-ray microscopes includes the VersaXRM-500, MicroXCT-200 and MicroXCT-400 for 3D X-ray computed tomography. The UltraXRM family of 3D X-ray microscopes includes the UltraXRM-L200 and UltraXRM-S for laboratory and synchrotron facilities, respectively.